

5R

QUICKIE®



Highly adjustable to elevate your mobility. The Quickie® 5R is designed to deliver superior rigid wheelchair performance day after day. As you and your environment evolve, or as you refine your riding capabilities, the 5R's frame and components adjust to elevate your mobility.

Offers two frame options. The longer Classic frame (C) and the shortened version (P). **Image 1.**

Swept-Back, Highly Adjustable Axle Plate Design. **Image 2.**

For an incredibly versatile backrest JAY J3 back with the back shell painted to match the chair. **Image 3.**

XTR Suspension Frame, Mountain wheels & Knobby tires. **Image 4.**



Image 1




Image 2



Image 3



Image 4




Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat Width	12" to 20"
Seat Depth	12" to 20"
Product Weight	10 kg / 22 lbs
User Weight Capacity	120 kg / 265 lbs
Frame angle	75°, 80° and 85°
Gravity center	Adjustable 1.5" - 4"
Accessories	Optional
Settings	In measurements
Rear suspension	Optional XTR

Available colours:



BENEFITS

For active people mainly spinal and motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

