

ACCESS

QUICKIE®



Its low seat-to-floor height allows access under tables and desks. It also allows the client to foot propel. If the client is using their hands to propel, the ACCESS tilt provides greater accessibility to the wheels while also promoting a stable and functional posture. Up to 25° of tilt and 55° of recline offers clinically effective pressure relief.

The ACCESS can be easily lifted & Portable. **Image 1.**

Rigid and Durable Frame. **Image 2.**

Transit Approved. **Image 3.**

Style in Motion. **Image 4.**



Image 1




Image 2



Image 3




Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

ITEM	EIR14
Width range	14" to 22"
Depth	14" to 22"
Product weight	From 14 kg Frame
Weight capacity	136 kg / 300 lbs
Tilt	-5° to 25°
Backrest angle	-30° to 50°
Backrest height	15° to 24"
Available colours:	



BENEFITS

Allows better access while allowing the user to self-propel.

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

