

# Advance Lift



 Joerns

The Hoyer® Advance® is compact, yet able to perform an outstanding range of transfers. Whether from the floor or onto a healthcare bed, the Advance copes with ease. The Advance's unique tapered leg design allows the lift to get even closer to the widest of obstacles. This is especially useful when encountering large chairs, bulky commodes, and wheelchairs. The tapered design also ensures that the resident feels safe and secure during the transfer.

Compact triangular folding design for easy transport and storage. **Image 1.**

Ergonomic design considerations and interaction points. **Image 2.**

Integrated safety system to help protect patients. **Image 3.**

Angled leg design ensures excellent access around furniture. **Image 4.**



Image 1



Image 2



Image 3



Image 4

# Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

## SPECIFICATIONS

Closed base	67 cm / 26.3"
open base	108 cm / 42.5"
Maximum height	186 cm / 73"
Minimum height	136 cm / 53.5"
base height	3 cm / 1.1"
Wheel size	3"
Weight capacity	155 kg / 341 lbs
Product weight	31.7 kg / 69.9 lbs
System	Electric (optional hydraulic)
Accessories	Transfer meshes



## BENEFITS

Essential in transfers for patients and caregivers in a comfortable and safe way.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

