

# Steel Wheelchair



Reliable and affordable, this steel wheelchair provides basic comfort and durability for both short- or long-term use. Ideal for: Patients in public or private settings looking for a classic, economic, and sturdy wheelchair.

Folding frame for easy transport and storage. **Image 1.**

Padded, folding armrests facilitate the user's lateral transfer. **Image 2.**

Optional 8" x 1" solid front casters. **Image 3.**

Optional elevating legrest. **Image 4.**

**ONE**




Image 1



Image 2



Image 3






Image 4



# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS	
Frame	Steel Folding
Seat width	16", 18" & 20"
Seat depth	16"
Product weight	18 kg
Weight capacity	120 kg / 265 lbs
Front wheel	8"x 2" (Optional 8" x 1")
Rear wheel	24" solid
Armrests	Folding
Foot support	Removable and adjustable Optional elevatable
Backrest	Standard fabric
HCPCS	K0001

	Footrest detachable and adjustable in height
	Rear wheels solid
	<b>Shipping Measurements</b> A 33" B 11" C 37" <b>Weight 21 kg.</b>

## BENEFITS

For people who require a light and safe transfer chair.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

