

Cushion Basic

JAY®



Gentle contour for basic needs with a soft, durable base gently contoured foam designed for the client at very low risk of skin breakdown requiring mild stability.

The soft Basic foam base is slightly contoured to provide support and increase forgiveness when. **Image 1.**

A moisture-resistant cover helps with incontinence management by protecting the basic foam base from moisture buildup. **Image 2.**

Water resistant materials easy to clean. **Image 3.**

Allows options in exact measurements in inches in width and depth. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Wheelchair cushions designed to provide the user with support, comfort and protection.

SPECIFICATIONS

Width range	14" to 24"
Depth	16" to 20"
Product weight	0.8 kg / 1.7 lbs (16"x16")
Height	5 cm / 2"
Weight capacity	14" - 20" 113 kg / 250 lbs 22" - 24" 136 kg / 300 lbs
Foam base	Pre-tuned foam
Outer cover	Anti-slip
Settings	In measurements

BENEFITS

Designed for users at low risk of skin degradation requiring mild stability

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

