

E-Pacer Transfer & Mobility



Its functionality and combination of lifting and gait training with large or highly dependent clients.

It can be raised high enough to support a 6'5" walking customer and then retracted to 32¾" for compact storage.

Image 1.

Turning the leg expander bar extends the frame of the E-Pacer to accommodate large wheelchairs.

Image 2.

The scale pairs with a tablet or phone via Bluetooth, the Gait Tracker app can track your client's weight-bearing progress.

Image 3.

The switch post attachment allows a caregiver or user to operate the electric lift without bending.

Image 4.



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
User measurement	77 cm to 193 cm / 30" to 76"
Trunk width	56 cm to 152 cm / 22" to 60"
Max user weight	160 kg / 350 lbs
Product weight	32.5 kg / 72 lbs
Use	Transfers and mobility
Accessories	Optional
Settings	No tools required
APP	Bluetooth, Gait Tracker app
Scale	Optional (factory)
Battery / Charger	Acid gel 24 V, 2,9 Ah / 100–240 V CA
Frame colors	

Required user measurement (Total height)	
	Odometer option
Twist release Twist lock Variable drag	
	Foot brakes and steering lock

BENEFITS

Excellent help for patients with different transfer needs and limited mobility.

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

