

# EIO Push Chair



Designed for children from early intervention to pre-teen ages, easy to fold and carry. The calf panel and backrest position your child correctly.

Legrest angle adjustment from +15 degrees to -55 degrees. **Image 1.**

Padded seat sides and backrest. **Image 2.**

Tomato EIO special rain cover. **Image 3.**

Folds compactly for transport - no need to remove rear wheels to fold. **Image 4.**



Image 1



Image 2



Image 3








Image 4



# Keeping you in mind. More options.

Clinically sound, innovative design and exceptional value.

SPECIFICATIONS	
Overall user height	76.5 cm - 127 cm / 30" - 50"
Seat width	14" / 35.5 cm
Seat depth	12" / 30.5 cm
Product weight	10 kg / 22 lbs
Backrest height	24" / 61 cm
Weight capacity	40.8 kg / 90 lbs
Recline angle	25° to 75°
Harness	Comfort Wrap 5 points
Head support	Including
Accessories	Optional
Color	Black base

Required user measurement (Total height)	
	Front wheels solid and removable 9 inches.
Rear wheels solid and removable 11 inches.	
	Compatible for Special Tomato seats

## BENEFITS

A special needs stroller that provides mild to moderate positioning support

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

