

# Early activity system

## LECKEY



It is designed to be used as part of an early intervention program to maximize a child's future potential.

1. Floor Mat
2. Wedge
3. Headrest
4. Sitting support
5. Positioning straps (pair)
6. Small support roll (pair)
7. Medium Support Rolls (Pair)



# Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

## SPECIFICATIONS

Age	Maximum 48 months
User weight	Maximum 15 kg
User height	Maximum 100 cm.
Floor mat	76 cm. Width 95 cm. Long
Sitting support	19 cm. Height 30 cm. Deep
Seat width	20 cm. minimum 25 cm. maximum
Small roll holder	6 cm. Width 60 cm. Long.
Medium Roll Support	8 cm. Diameter 70 cm. Long.
Trunk Wedge	27 cm. Width 20.8 cm. Long. at 7.5°
Straps pos.	8 cm. Width 65 cm. Long.



## BENEFITS

Suitable for babies and toddlers from birth to 4 years.

## Discover our products **ONE**



### DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

