

Fit Backrest



 **JAY**®

The Jay Fit backrest provides multiple opportunities for positioning and growth. With its anatomical shape and multiple settings, the Jay Fit backrest gives the patient greater comfort and independence, providing trunk stability and proper positioning. Thoracic laterals and head supports contribute to the correct positioning of the user.

Allows multiple options for growth and change over time. **Image 1.**

Now Available with Whitmyer® Headrests and JAY Positioning Supports. **Image 2.**

New Swing-Away Lateral Design. **Image 3.**

Mounting Hardware options. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Wheelchair cushions designed to provide the user with support, comfort and protection.

SPECIFICATIONS

Width range	10" a 18"
Height	10" a 22"
Depth	8" a 20"
Product weight	2.5 kg / 5.5 lbs
Weight capacity	113 kg / 250 lbs
Adjustment angle	Allowed
Material	Aluminum and foam
Cover	Microclimatic 3DX
Settings	In measurements
Accessories	Optional

Configurations with
JAY cushions



BENEFITS

For active users with lack of sensitivity and stability

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

