

# Frontier V6

**magic**  
Mobility



All-terrain mid-wheel drive power wheelchair provides tremendous off-road and curb-climbing capabilities. Its low pressure tires provide a smoother ride and give the Frontier V6 increased climbing ability. The Frontier V6 is the ideal outdoor wheelchair, yet its size makes it a capable indoor performer, too!

With a 12" seat lift, the C-Me electric lift function. (5) Travel at 5 km/h with 12" seat elevation. **Image 1.**

Powerful LED lights Installed on rubber mounts, these lights are flexible and impact resistant. **Image 2.**

The Frontier V6 Hybrid Wheel and Tire Package. **Image 3.**

Offers a wide range of options and accessories. **Image 4.**



Image 1




Image 2



Image 3



Image 4




# Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

## SPECIFICATIONS

Seat width	13" up to 23"
Seat depth	13" up to 24"
Product weight	From 145 kg
Weight capacity	154 kg
Backrest height	16.5" up to 24"
Maximum range	33 km
Maximum speed	8 kmh
Wheels	8" and 14"
Turning radius	58 cm (without footrest)
Total width	27.3" and 28"



- Seat lift and tilt 12" / 45°
- 95-180° Legrest elevation
- Tilt 50°
- 95-160° Recline

## BENEFITS

For users with partial or total reduced mobility of lower limbs.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

