

HTS

Hygiene and Toileting System



The HTS is yet another fantastic innovation. The creative and highly adjustable design gives children comfort, security, while offering caregivers ease and simplicity. Simple adjustabilities and ease of cleaning. The quality of this product is by far the best on the market.

The Tri-Frame is a new design that better In any toilet model, round or elongated front, it allows the user to use the same toilet as the rest of the family. *Image 1.*

Various models of bases, fixed and mobile (tilting in space with stirrup). *Image 2.*

Three models with seat and back padding. *Image 3.*

Large wheels improve mobility when overcoming higher thresholds or other obstacles. They add approximately 1" to the height from the floor. *Image 4.*



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. **More options.**

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
Total user measurement	76 cm - 188 cm / 30" - 74"
Frame measurements	Small, Medium y Large
Max user weight	34 kg - 113 kg / 75 lbs - 250 lbs
Product weight	4.6 kg - 14.5 kg / 10.2 lbs - 32 lbs
Use	WC, shower chair and mobile base
Accessories	Optinal
Settings	No tools required



BENEFITS

Suitable for patients with neuronal and psychomotor damage or cognitive psychic disability.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

