

HiLite XTR

RGK



Made from the highest quality materials and hand crafted with care. The HiLite XTR is functional, lightweight and comfortable. The HiLite XTR is the perfect Titanium chair to experience made to measure for the first time, whether a new user or an experienced one. With over 30 measurements to be tailored to the user, the clinical benefits and comfort levels are life changing. A lightweight Titanium grade 2 frame and a huge array of options to choose from, ensures superior functional performance in this classic design wheelchair.

The wheelchair frame is made from commercially pure Titanium grade 2.

Image 1.

With over 30 measurements, can be perfectly tailored to suit the user's disability and clinical requirements.

Image 2.

Fixed or folding backrest, fixed height or adjustable height backrest.

Image 3.

Airtech upholstery system. **Image 4.**



Image 1



Image 2



Image 3




Image 4

Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame Material	Titanium grade 2
Frame style	Rigid
Seat Width	10" - 18" / 25 cm - 46 cm
Seat Depth	10" - 18" / 25 cm - 46 cm
User Weight Capacity	125 kg / 275.6 lbs
Product Weight	From 7.5 kg / 16.5 lbs
Folding system	Q-lock folding system
Available colours	25 options



Rear wheels
Quick release

BENEFITS

For active people mainly motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

