

# Hippocampe All-terrain Chair



Since 2002, the Hippocampe wheelchair offers individuals with disabilities accessibility to difficult terrains, impossible in a standard wheelchair: beach, sand, swimming area, hiking trails, snow, mountain.



# Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

## SPECIFICATIONS

Frame	Stainless steel and anodized aluminum
Wheels	Optional for pool or beach
Seat width	39.4 cm /15.5"
Seat depth	36.8 cm /14.5"
Backrest height	36.8 cm /14.5"
Total depth	72.4 cm /28.5"
Weight capacity	130 kg / 286 lbs
Accessories	Optional
Frame colors	Blue and pink



## BENEFITS

Excellent help for patients with different transfer needs and limited mobility.

## Discover our products



### DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

