

Horizon Stander



LECKEY

A revolutionary 3-in-1 stander, it is designed for users from kids to adults with complex postural needs. This electric-operated standing system provides outstanding postural support and pelvic stability in either prone, supine, or upright positions

Include chest, hip, knee, and footplate supports to provide easy adjustability and optimal standing positioning.

Image 1.

The Leckey Horizon Stander is three standers built into one. Usable as an upright, supine, or prone stander.

Image 2.

The versatile and easy-to-use angle adjustment a standing position that has shown to have many benefits to overall health and everyday life.

Image 3.

The Leckey Horizon Stander is perfect for the classroom and home settings, allowing eye-to-eye contact and freedom for hands to move and develop motor skills. **Image 4.**



Image 1



Image 2



Image 3




Image 4

Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
Sizes	3
Age	4 to 14 years
User height	100 to 180 cm / 39.4" to 71"
Max user weight	100 kg / 220 lbs
Angle adjustments	0° -90° - Prone 0° -90° - Supine
tray angle	0° -40° - Prone +/- 20° / 30° - Supine
Fabric options	Grey, Blue and Pink

Required user measurement (total height)



BENEFITS

Robust and durable 3-in-1 stabilizer for users from 4 years old to adults.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

