

J3 Backrest

Aluminum and carbon



The most versatile JAY back for all shapes, sizes, and lifestyles. The JAY® J3™ backrest product line offers a multitude of width, height, and contour depth back shells to fit the back to the client. Each back is lightweight and adjustable, with numerous hardware choices allowing for easy installation.

The J3 Plus offers pelvic and trunk stability for the bariatric client.

Image 1.

More sizes. More styles. PA, PL, PD, PDL. **Image 2.**

The most popular 6" contoured backrest, PD, is now available in over 40 sizes. **Image 3.**

Offers a carbon fiber shell option for the ultimate in lightweight performance. With such a light weight from 2.2 lbs. **Image 4.**



Image 1




Image 2



Image 3



Image 4



Keeping you in mind. More options.

Wheelchair backrest designed to provide the user with support, comfort and protection.

| SPECIFICATIONS | |
|------------------|--|
| Width range | 12" to 26" |
| Height | 6.5" to 24" |
| Depth | 2.2" to 6" |
| Product weight | 0.8 kg / 1.7 lbs (without supports) |
| Weight capacity | 12"-20" 136 kg / 300 lbs 22"-26" 226 kg / 500 lbs |
| Adjustment angle | Allowed |
| Material | Soft foam |
| Sheath | Microclimatic 3DX |
| Base | Aluminum or Carbon |
| Settings | Height, width and depth |
| Accessories | Supports and harness |



BENEFITS

Backrests designed for wheelchairs for perfect support

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

