

Compact Scooter

limosa 



 **ONE**

Compact enough to travel with you but rugged enough to take outdoors, Limosa® scooter combines the performance of a sport scooter with the portability of a travel scooter.

Front and rear suspension. **Image 1.**

Detachable 5-piece design fits in most trunks or rear seats. **Image 2.**

Height adjustable seat and rudder adjustable ensure you will have a comfortable ride. **Image 3.**

The external charger makes charging easy. Directly connect or carry box of the battery to a place to charge. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

| | |
|----------------------|---------------------|
| Frame | Detachable 5 pieces |
| Seat Width | 17" |
| Seat Depth | 15" |
| Backrest height | 15" |
| Product weight | 49 kg / 108 lbs |
| Weight capacity | 136 kg / 300 lbs |
| Wheels | 8"x 2.5" Solid |
| Base to floor height | 4" |
| Armrests | Folding and padded |
| Speed | 6 mph |
| Color | Green oil |



Shipping dimensions

A 47"

B 20"

C 21"

Weight 55 kg.

BENEFITS

Designed to be driven on sidewalks and pedestrian areas.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

