

# MPS Car seat



The special seat offers a solution if you plan to use in the car or at home. The seat has been designed to offer comfortable positioning and proper support in the car or in the optional base.

The Seat has anti-microbial properties, is impermeable to fluids, as well as latex-free and peel- and tear-resistant.

**Image 1.**

30° Hip-Flex® Adjustment for optimal pelvic positioning. **Image 2.**

Seat shell in two sizes, small and large. **Image 3.**

Shoulder Strap Cover. **Image 4.**



Image 1




Image 2



Image 3








Image 4



# Keeping you in mind. More options.

Clinically sound, innovative design and exceptional value.

SPECIFICATIONS	
User height	Size 1 to 4 75 - 185cm / 28" - 62"
User weight	Size 1 to 4 36 up to 59 kg / 80 - 130 lbs
User age	From 9 months to adults
Seat width	7"-13"
Deep seat	6" - 17"
Harness	5 points car seat grade
Accessories	Optional
Colors	Aqua, Cherry, Lilac and Gray 

Required user measurement (Total height)	
	Safety Attachment Straps For use outside of a vehicle
Tray with Risers	
	Mobile Floor Base

## BENEFITS

They are intended for children with mild to moderate positioning needs.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

