

MaxLite

RGK



Only the best materials are used in your MaxLite. Titanium is famous for having all of the positives of aluminium with the additional bonus of a natural shock absorption by the material itself. The utmost best performance of your chair is ensured by only using elements produced by market leaders, alongside a staggering 19 quality checks throughout the build, from measure to handover.

The wheelchair frame is made from commercially pure Titanium. **Image 1.**

With over 30 measurements, can be perfectly tailored to suit the user's disability and clinical requirements. **Image 2.**

Choose an accent colour to suit your style. **Image 3.**

Airtech upholstery system. **Image 4.**



Image 1



Image 2



Image 3



Image 4

Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame Material	Titanium
Frame style	Rigid
Seat Width	10" - 17" / 25 cm - 44 cm
Seat Depth	10" - 17" / 25 cm - 44 cm
User Weight Capacity	100 kg / 220 lbs
Product Weight	From 7.5 kg / 16.5 lbs
Accessories	Optional
Available colours	25 options



Rear wheels
Quick release

BENEFITS

For active people mainly motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

