

# NITRUM

**QUICKIE®**



For active users looking for different driving performance, innovative details and a wide variety of options, Nitrum creates the perfect combination of fit and performance.

The Nitrum Hybrid offers a reinforced front frame without compromising its lightweight construction. **Image 1.**

For better positioning and functionality of the chair, it can be configured with JAY backrests and cushions. **Image 2.**

Made of 7020 oval aluminum with ShapeLoc technology and removable frame. **Image 3.**

The first integrated LEDs with removable battery. **Image 4.**



Image 1




Image 2



Image 3



Image 4



# Pensando en ti. Más opciones.

For those who require cutting-edge technology to keep up with the demands of their active life.

## SPECIFICATIONS

ITEM	EIR12
Seat width	12" a 20"
Seat depth	12" a 20"
Product weight	From 5 kg the frame
Weight capacity	124 kg / 275 lbs
Camber	0°, 3° y 6°
Gravity center	Adjustable 0" up to 4"
Accessories	Optional
Available colours:	

A unique twist for folding the backrest Twist-Lock



Reinforced frame

Patented adjustable wheel system



Custom colors included

## BENEFITS

For active people mainly spinal and motor injuries.

## Conoce nuestros productos



### DISEÑADOS PARA NUESTROS TERRITORIOS

Nuestros productos están diseñados para satisfacer las necesidades de cada uno de nuestros territorios, ofreciendo calidad y precios accesibles.



### VERSÁTILES Y FUNCIONALES

Contamos con diversidad de productos que van desde una gran gama de sillas de ruedas manuales, motorizadas y scooter de movilidad.

Para más información sobre cómo obtener este producto  
**Contacta a tu representante local de Loh Medical**

Encuentre su representante local en [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.115

