

DYNAMIC PACER Gait Trainer

 rifton.



The Pacer gait trainer is the solution for those who are starting to walk and need assistance in the process of rehabilitation and strengthening of muscles and bones safely.

The standard top frame is a great option for customers who do not need dynamic movement. **Image 1.**

The dynamic upper frame allows movement to provide dynamic weight shifting and dynamic body weight support to encourage natural gait patterns. **Image 2.**

The standard base is designed for use both indoors and outdoors. **Image 3.**

Utility base for navigating over doorways, lawns, gravel paths or chipped play areas. **Image 4.**



Image 1



Image 2



Image 3







Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
Elbow to floor meas.	39 cm - 127 cm / 15.5" - 49"
Frame measurements	Mini, Small, Medium, Large y X Large
Max user weight	34 - 113 kg / 75 - 250 lbs
Product weight	4.6 - 14.5 kg / 10.2 - 32 lbs
Use	Indoor and outdoor
Accessories	Optional
Settings	No tools required
Frame colors	

Required user measurement (Elbow to the floor)	
	Odometer option
Twist release Twist lock Variable drag	
	Foot brakes and steering lock

BENEFITS

For users who need to improve their ambulatory skills and develop weight support

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

