

PAL Seating System

LECKEY



Reinforced chair that is easy to use at school, home or rehabilitation centers. For mild or moderate postural needs, it improves sitting stability, which allows attention periods to be prolonged and facilitate manipulation activities and fine motor activity.

Two-point lap belt and side pads for lateral stability. **Image 1.**

Height adjustable armrests, legs and footplate. **Image 2.**

Base and footrest options available, including castors, rockers and stabilisers. **Image 3.**

Adjustments Made Easier Than A-B-C. **Image 4.**



Image 1



Image 2



Image 3



Image 4




Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.


SPECIFICATIONS

Sizes	4
Aproximated age	1 to 12 years
Max user weight	30 to 50 kg /66 up to 110 lbs
User height	75 to 160 cm / 29.5" up to 63"
Seat width	15 up to 38 cm / 6" up to 15"
deep seat	15 up to 40 cm / 15" up to 15.7"
Backrest height	26.5 up to 45 cm / 10" up to 17.7"
Backrest recline	From -95° to +105°

Colors



Required user measurement (total height and knee to heel)



BENEFITS

Pediatric postural seating system, with mild to moderate postural needs.

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

