

PRONESTANDER



LECKEY

In addition to prone angle adjustment, the stander covers a wide range of growth and has adjustable chest, hip and knee supports, to ensure it meets a wide range of clinical needs.

Easy to set up and adjust for multiple users. **Image 1.**

Five incremental secure prone positions possible, from 85° to 45°. **Image 2.**

Supplied as standard with individually adjustable chest, hip, knee and sandal supports. **Image 3.**

Enables de-rotation of the hips as well as alignment of the pelvis and centre of gravity over the base of support. **Image 4.**



Image 1



Image 2



Image 3




Image 4


Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Sizes	3
Age	From 1 to 18 years
User height	From 80 to 185 cm / 31.5" to 72.8"
Max user weight	90 kg / 198 lbs
Prone angles	85°, 75°, 65°, 55°, 45°
Fabric options	Size 1 Prairie Green Size 2 and 3 Blue Mondeo

Required user
measurement
(total height)



BENEFITS

Pediatric postural seating system, with mild to moderate postural needs.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

