

Q700-UP M

 QUICKIE®



Providing the highest performance and smartest seating technology of all our power wheelchairs. Whether you want to drive off-road, at home, in the shops, lie down or stand up, the Q700-UP M lets you do all this and much more.

It offers the recommended option for JAY cushions and backrests tailored to the user's specific needs. **Image 1.**

Switch-It Specialty Controls and 6 Memory Positions. **Image 2.**

Allows you to stand whenever you want and offers many benefits for your health and lifestyle. **Image 3.**

Allows for indoor and outdoor use thanks to the base accessories, front and rear wheels. **Image 4.**




Image 1




Image 2




Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	16" to 20"
Seat depth	17" to 22"
Backrest height	22" to 28"
Weight capacity	120 kg / 265 lbs
Product weight	From 198 kg >
Wheel system	Spyder Track 14" and 6"
Features to choose	Standing, Tilting, Lifting Recline, Elevation of feet
Motors	4 poles
Speed	9.6 kmh standard
Bluetooth	Optional

Available colours:



BENEFITS

For active people mainly spinal and motor injuries.

Electronic functions that can be configured and customized



2 sealed batteries

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

