

R90

Empulse



Travel longer distances without getting tired, climb ramps with ease, and reduce shoulder strain. The sleek, simple, compact design allows the R90 to fit under your chair. It's discreet and intuitive, and maintains the size and turning radius of your chair.

A push of a button elevates the wheel for zero rolling resistance. **Image 1.**

Start, stop, and control speed with the push of a button. **Image 2.**

Sunrise Intelligence App lets you select performance profile, top speed, and displays battery charge. **Image 3.**

V-Lock attachment system secures the R90 to most folding and rigid wheelchairs. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

ITEM	EISP8
Maximum range	Standard Battery (2.5Ah): 9.8 miles Extended Battery (5Ah): 19.5 miles
Maximum speed	5.5 mph (adjustable)
Battery	36V Lithium-Ion, 2.5Ah/5Ah
Weight capacity	275 lbs
Product weight	R90 without battery: 13.0 lbs Standard Battery: 1.5 lbs Extended Battery: 2.6 lbs



BENEFITS

Powered Assist Device for Folding & Rigid Frame Wheelchairs

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

