

# Ranger 2 Pool Lift

  
Aqua Creek Products



The Ranger 2 upholds the same standards for durability and reliability as all Aqua Creek pool lifts. It is designed to be independently operated by the user, allowing full control of the lift with an easy to use 2 button handset. The Ranger 2 is ideal for hotel and motel pools and also a great choice for backyard pool enthusiasts.

Standard anchor. **Image 1.**

Standard anchor with quick connect knobs. **Image 2.**

Pavement anchor. **Image 3.**

Wood Deck Anchor. **Image 4.**



Image 1



Image 2



Image 3



Image 4



# Keeping you in mind. More options.

The leaders in recreation, fitness and accessibility, safe and guaranteed.

## SPECIFICATIONS

Frame Finished stainless steel  
in powder paint

Seat width 45.7 cm / 18"

Seat depth 43.2 cm / 17"

Backrest height 44.3 cm / 18"

Arm support Folding

Foot supports Adjustable and belt

Weight capacity 160 kg / 350 lbs

System Electric

Accessories Optional

Colors



Solar charger  
helps extend  
battery life



Outdoor cover

Match your branding  
or decoration with  
custom color



## BENEFITS

Manufactured for pools and spas, they are ideal for both commercial and residential use.

## Discover our products **ONE**



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product  
Contact your local Loh Medical representative

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

