

Squiggles Seating system



LECKEY

This pediatric seating and positioning system provides unparalleled postural positioning of the pelvis, torso, head, and legs while giving the child the greatest possible freedom of movement.

The Hi-Low base allows for the child's height to be easily changed via a gas compression foot pedal for eye-to-eye interaction during work, rest, and play.

Image 1.

Large variety of adjustment possibilities, the Squiggles Seating System is able to grow with your child.

Image 2.

Pelvic stability is the basis for comfort and optimal postural support in every therapy chair.

Image 3.

Leckey Seating on ZIPPIE Mobility Bases.

Image 4.



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.


SPECIFICATIONS

Chassis base	Hi-Low
Aproximated age	From 1 to 5 years
User height	From 75 to 110 cm / 29.5" to 43.3"
Max user weight	22 kg / 48.4 lbs
Backrest recline	From -10 to 25°
Seat tilt	From -10 to 30°
Footrest angle	From -10 to 10°
Postural unit weight	6 kg / 13 lbs
Hi-Low base weight	9 kg / 20 lbs

Colors



Required user
measurement
(total height)



BENEFITS

Designed for children and young adults who need optimal positioning and pelvic stability

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

