

SR45



The simple solution for rotation-in-space technology. The QUICKIE® SR45™ tilt-in-space wheelchair features 45° of Intelligent Rotation in Space Technology and a simplified mix of options to provide proper positioning and comfort at an outstanding value.

45° of intelligent rotation in space technology. **Image 1.**

Durable frame. **Image 2.**

Power tilt, attendant control, toggle switch. **Image 3.**

Quickie's top options. **Image 4.**



Image 1



Image 2



Image 3




Image 4




Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Frame	Steel
Seat Width	14" to 20"
Seat Depth	14" to 22"
Backrest height	15" to 24"
User Weight capacity	120 kg / 265 lbs
Product weight	18 kg / 39 lbs
Tilt degrees	Up to 45°
Tilt system	Tilt in space
Accessories	Optional
Settings	In measurements

Available colours:



Pedal system
to activate
tilt and brake



Intuitive
foot support designs



BENEFITS

A special needs chair that provides tilting and positioning support.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

