


Suitcase Ramp



This portable ramp, with its single-fold design, sets up quickly and is ideal for scooters and wheelchairs. Made of lightweight aluminum with an easy-to-carry design, the SUITCASE® Singlefold Ramp is fitted with ergonomically-designed handles and a full-length hinge for additional strength and stability. Designed for high-traffic use, this durable ramp is available in 2' to 8' lengths and features an extruded slip-resistant tread for superior traction.



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Material	Aluminum
Broad	75 cm / 29.5"
Deep	61 - 243 cm / 2' - 8' Ft.
Product weight	3.6 - 13 kg / 12 - 30 lbs
Weight capacity	362 kg / 800 lbs
Surface	Anti-slip
Creased	In two
Extension	Optional



BENEFITS

Suitable for home, vehicle or travel and offer temporary or long-term solutions

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

