

# Scout Excel



The Scout Excel pool lift features an adjustable height seat pole, upgraded electronics and a more robust base. These upgrades allow for a weight capacity of 375 pounds. This pool lift will accommodate most existing anchors by utilizing adapter sleeves. The optional pull-out leg rest may be used for greater clearance on raised pool and spa applications.

Standard anchor. **Image 1.**

Semi-recessed anchoring. **Image 2.**

Pavement anchor. **Image 3.**

Wood Deck Anchor. **Image 4.**



Image 1



Image 2



Image 3





Image 4



# Keeping you in mind. More options.

The leaders in recreation, fitness and accessibility, safe and guaranteed.

SPECIFICATIONS	
Frame	Finished stainless steel in powder paint
Seat width	45.7 cm / 18"
Seat depth	45.7 cm / 18"
Backrest height	45.7 cm / 18"
Arm support	Folding
foot supports	Adjustable and belt
Weight capacity	170 kg / 375 lbs
Rotation	360°
System	Electric
Accessories	Optional
Colors	

Solar charger helps extend battery life	
	Outdoor cover
Match your branding or decoration with custom color	

## BENEFITS

Manufactured for pools and spas, they are ideal for both commercial and residential use.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product**  
**Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

