

# Soft-Touch Liners



Special  
**tomato**

They offer additional support and comfort and are recommended for children with mild to moderate sitting support needs. It is made of a soft-touch material that is latex-free and antimicrobial. The seat and back cushions can be used together or separately.

Soft contoured seat cushion and/or back cushion. **Image 1.**

They can be attached to most seating systems. **Image 2.**

It has been designed and molded with a "hinge" to move comfortably and support the body while reclining. **Image 3.**

They are easily transported. **Image 4.**

[www.lohmedical.com](http://www.lohmedical.com)




Image 1



Image 2




Image 3



Image 4



# Keeping you in mind. More options.

Clinically sound, innovative design and exceptional value.

## SPECIFICATIONS

Sizes	Size 1 to 5
Seat cushion width	20 - 40.6 cm / 8 - 16"
Deep seat cushion	20 - 40.6 cm / 8 - 16"
Backrest width	18.75 - 35.5 cm / 7.5 - 14"
High backrest	25 - 40.6 cm / 10 - 16"
User weight	9 up to 91 kg / 20 - 200 lbs
Colors	 Cherry  Dark gray  Lilac  Aqua

3 Point Pelvic Cushion Harness and Trunk Harness.



Built-in abductor for seat cushion

For special recliner support system



## BENEFITS

They are intended for children with mild to moderate positioning needs.

## Discover our products



### DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

