

Splashy Big Multipurpose Bath Seat



Portable bath seat that makes showering easier and much more fun for families of children with special needs. It's a portable, lightweight, supportive seat that works for bathing or messy play wherever you go.

Splashy is not just for the bath or shower, it can be used at the beach, in a baby pool, or for play activities. **Image 1.**

Adapts to your child's postural support needs, with the option of 5 or 3 point pelvic support. **Image 2.**

4 removable positioning supports for use anywhere on the seat where support is required. **Image 3.**

26 recline options, from upright at 106 degrees to fully reclined at 140 degrees. The Floorsitter can also be used with your GoTo Seat. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Backrest height	72 cm / 28.5"
Seat width	45 cm / 17.5"
Deep seat	40 cm / 15.5"
Backrest width	45 cm / 17.5"
Max user weight	55 kg / 121 lbs
Age range	8 to 16 years
Colors	Seat and backrest Blue Accessories, pink and green Gray base

The Floorsitter can also be used with your GoTo Seat.



BENEFITS

Suitable for patients with neuronal and psychomotor damage or cognitive psychic disability.

Discover our products



DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

