

Strive



The Strive is a lightweight, easy-fold stroller with positioning options. It is designed for travel and provides upper body support and full body positioning. Optional features, adjustments and accessories meet the needs of children with disabilities.

Folding frame easy to carry and store. **Image 1.**

Accommodates growth with its adjustable seat depth. **Image 2.**

Five-point harness and pelvic belt. **Image 3.**

Wheel Lock and traffic system. **Image 4.**




Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
Seat width	12", 14", 16" and 18"
Seat depth	11.5" up to 17"
seat angle	80°, 90° and 95°
Tilting	20°
Backrest height	24"
Wheels	8" front, 10" rear
Backrest	Adjustable tension
Product Weight	15 - 16 kg / 33 - 36 lbs
Weight capacity	34 - 113 kg / 75 - 250 lbs
Optional	Headrest, Canopy, Basket under the seat and anti-tip
Fabric colors	Red and blue



BENEFITS

Available in four sizes with adjustable seat depth for years of growth!

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

