

TS



ZIPPIE's original pediatric tilt-in-space wheelchair. The ZIPPIE® TS™ is the original, premier pediatric tilt-in-space wheelchair. Available with either a folding or rigid frame, built-in growth, and compatible with many seating and positioning solutions, the TS™ has the widest variety of options in the tilt-in-space market.

45° Traditional Tilt Technology. **Image 1.**

Folding and Rigid Frame Options. **Image 2.**

Growth Capabilities. Up to 2" of seat width and up to 5" of seat depth adjustments. **Image 3.**

ZIPPIE TS with Leckey Seating. **Image 4.**



Image 1



Image 2



Image 3



Image 4




Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Frame	Aluminum
Frame type	Foldable / Rigid
Seat width	10" to 18"
Seat depth	13" to 20"
User Weight capacity	75 kg / 165 lbs
Product weight	13.5 kg / 29 lbs
Tilt	Up to 45°
Accessories	Optional
Settings	In measurements

Available colours



BENEFITS

A special needs chair that provides tilting and positioning support.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

