

Tiga Jnr

RGK



The Tiga Junior is specifically designed to your child's individual measurements, and can be personalised through a choice of colours and embroidery. The unique intelligent growth system provides maximum flexibility and adjustability to grow with your child. Give your child the best start in life.

Intelligent growth system with adjustable features. **Image 1.**

Airtech Technology. Range of ergonomic backrest options. **Image 2.**

Low ride castor forks for optimum performance. **Image 3.**

A full range of active & assist options. **Image 4.**



Image 1




Image 2



Image 3




Image 4

Keeping you in mind. More options.

Quality products with the required accessories

SPECIFICATIONS

Frame Material	7020 Aluminum
Frame Style	Rigid
Weight Capacity	70 kg / 154 lbs
Product Weight	From 7.5 kg / 16.5 lbs
Seat Depth	Maximum 15.7" / 40 cm
Seat Width	Maximum 14" / 36 cm
Growth System	Intelligent growth system

Add personal embroidery on the backing.



Rear wheels
Quick release

BENEFITS

For active people mainly motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

