

WeGo



Lightweight, portable and stylish like any typical stroller and you can also place your GoTo in it, so your little one gets all the support and comfort for strolling.

It's really easy to get your GoTo Seat in and out of the WeGo. **Image 1.**

It reclines to 140° and they'll still have the support of the GoTo Seat. **Image 2.**

Special rain cover. **Image 3.**

Folds compactly for transport - no need to remove rear wheels to fold. **Image 4.**




Image 1



Image 2



Image 3



Image 4




Keeping you in mind. More options.

Clinically sound, innovative design and exceptional value.

SPECIFICATIONS

Total user height	116 cm / 45.6"
Seat width	35 cm / 13.5"
Deep seat	26 cm / 10"
Max user weight	30 kg / 66 lbs
Product weight	11 kg / 24.2 lbs
Age range	1 to 8 years
Accessories	Optional
	Go To Seat Compatible

Required user
measurement
(Total height)



BENEFITS

A special needs stroller that provides mild to moderate positioning support

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

