

X'CAPE



It is designed with the goal of helping children with different mobility needs move to the next level of independence. One of the lightest pediatric folding wheelchairs at just 19 pounds, it allows your child to spend less energy on propulsion.

Pediatric specific design, 2" rise in the crossbar and 3" in the seat rail. **Image 1.**

Frame with XLOCK folding spreader for performance and comfort. **Image 2.**

Z-finity® lower extremity positioning system. **Image 3.**

VersaRail™ for unlimited seating and accessory options. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS	
ITEM	EIFZ1
Frame	Aluminum
Frame type	Foldable / Rigid
Seat width	8" to 16"
Seat depth	8" to 16"
Weight capacity	75 kg / 165 lbs
Accessories	Optional
Settings	In measurements
Available colours:	

Lightweight and foldable frame easy to transport and store



BENEFITS

Technology of a rigid frame and folding frame for active and growing children.

Discover our products



DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

