

Compact Power Chair

FELIX^e



 **ONE**

A modern power wheelchair with a unique front-folding design that makes storage and transport easier than ever. Sleek, cushioned, and ready for indoor and outdoor use. Ideal for Patients with active lifestyles who value compact solutions without compromising on performance.

Folding frame for easy portability and storage if the user has reduce space.

Image 1.

The easy-to-use 4-button PG joystick extremely sensitive and controls dual Linix engines. **Image 2.**

Linix Dual Motors. **Image 3.**

Padded and folding armrests. **Image 4.**



Image 1



Image 2



Image 3




Image 4




Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Aluminum, Folding
Seat Width	18"
Seat Depth	17"
Backrest height	18"
Backrest type	Padded fabric
Product weight	36 kg / 79.4 lbs
Weight capacity	136 kg / 300 lbs
Wheels	12" solid 8"x 2" solid
Motors	200W*2 LINIX
Speed	7 kmh
Battery	Lithium 20Ah, HP 24V charger
Color	Black
HCPCS	K0013



Shipping dimensions

A 28"

B 16"

C 35"

Weight 42 kg.

BENEFITS For users with reduced mobility partial or total lower limbs.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

