

Q100 R

QUICKIE®



Designed to challenge the standards of indoor and outdoor power wheelchairs, the Q100 R is ultra-compact, offering agility indoors and the stability and performance you'd expect on city streets. Smart design, stylish appearance.

Weighing from 95 kg, it is very easy to lift and transport. **Image 1.**

PG VR2 Joystick Controller, high compatibility and easy installation. **Image 2.**

Compact at 21" in total width, practical and easy to drive. **Image 3.**

New, smaller center-mount options reduce the turning radius by 4". **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	17" to 19"
Seat depth	16" to 20"
Backrest height	20"
Weight capacity	125 kg
Product weight	From 95 kg >
Wheel system	8" and 12.5"
Motors	2 poles
Speed	6 kmh
Turning radius	63"

Available colors:



2 sealed batteries
55 Amp



BENEFITS

For active people mainly spinal and motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2026.02.21E

