

Aluminum Active wheelchair



 ONE



Image 1



Image 2



Image 3



Image 4

This active aluminum wheelchair features a rigid frame with maximum responsiveness and a modern style. Sleek and lightweight, it offers excellent maneuverability for dynamic patients. Ideal for: Patients with active lifestyles seeking a sporty, compact chair with high performance.

Adjustable center of gravity. **Image 1.**

Includes removable side clothing protectors with mudguard, adjustable in height and depth. **Image 2.**

Compact side brakes. **Image 3.**

Anti-tip wheels. **Image 4.**

Other features:


- Black frame Nylon seat and backrest
- Fabric backrest with adjustable 4-band tension
- Manual reclining backrest from -5° to 15°
- Center of gravity adjustable in 3"
- Aluminum footrest plate with foldable depth and angle adjustment
- Adjustment on the forks in 2"
- Armrests (Optional)
- Adjustable seat height
- Inflator included

Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Aluminum, powder coating
Seat width	16" & 18"
Seat depth	16"
Backrest height	14" - 16"
Total width	63 cm
Total depth	83.5 cm
Front wheels	5"x 1.5" solid aluminum wheel
Rear wheels	24" pneumatic, quick release
Weight capacity	110 kg. / 242.5 lbs.
Product weight	14 kg. 9 kg. without rear wheels
Optional accessories	Tubular Arm Supports
AAA AAP	18" - 20" 14" - 17"



Shipping dimensions

A 33"

B 25"

C 24"

Weight 17 kg.

BENEFITS

For people who require a chair for daily use, in a light and safe way.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2026.04.11E

