

Reclining Chair



This reclining wheelchair offers full back support, including a padded headrest and elevating leg rests for enhanced positioning. Its foldable parts make transport and storage easier. Ideal for: Patients requiring full body support in a low-cost and transportable solution.

Folding frame for easy transport and storage, with a locking bar for backrest recline. **Image 1.**

Footrests that elevate from 50° to 0°, are height-adjustable, foldable, and removable. **Image 2.**

Backrest reclines from 85° to 175°. **Image 3.**

Quick-release rear wheels for easy transportation. **Image 4.**




Image 1




Image 2



Image 3



Image 4



Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame structure	Steel Silver vein
Seat Size	16"x16", 18"x16" & 20"x17"
Product weight	26 kg
Weight capacity	120 kg / 264.5 lbs
Front wheels	8"x2"
Rear wheels	24" pneumatic quick extraction
Armrest	Detachable
Footrest	Elevatable, foldable and removable
Backrest height	20" / With head support 31"
Head support	Including



Shipping dimensions

A 36"
B 12"
C 38"
Weight 30 kg.

BENEFITS

A special needs chair that provides recline and positioning.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2026.04.11E

