

Power Chair rigid frame

FELIX[®]



The Felix[®] R combines comfort, postural support, and precise control in a motorized wheelchair with a rigid frame for indoor and outdoor use. The joystick allows for easy handling, while its ergonomic backrest and headrest improve the user's positioning. It also features a seat that adjusts to body measurements, allowing for personalized and comfortable use throughout the day.

Manual recline from 90° to 130°.

Image 1.

Easy-to-use joystick. **Image 2.**

8.5" x 2.5" solid black front wheels.

Image 3.

13.5" x 2.5" solid black rear wheels.

Image 4.

ONE



Image 1



Image 2



Image 3




Image 4




Keeping you in mind. More options.

Quality products with the required accessories at a great price.

SPECIFICATIONS

Frame	Rigid Blue steel
Seat Width x Depth	18" x 17"
Backrest height	20.5"
Backrest type	Padded Fabric
Product weight	102.4 kg / 226 lbs
Weight capacity	136 kg / 300 lbs
Wheels	13.5" x 2.5" Solid 8.5" x 2.5" Solid
Motors	2 x 3402W, Electromagnetic Brake
Max. Speed	10 kmh
Endurance mileage	30 km
Battery	12 V, 38 Ah Lead-Acid
Suspension	Rear
HCPCS	K0013



BENEFITS For users with reduced mobility partial or total lower limbs.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2026.04.11E

